

Cuticular

REVERSE CUTTING

Ethicon	Medtronic (Covidien)	Needle Length (mm)
1/2 Circle		
M-2		13
1/4 Circle		
C-2	C-1*	13
C-3	C-1*	13
FS-2	C-13	19
FS-1	C-14	24
FSL	C-15	26
FSLX	C-16	30
LR	GS-18	75
1/2 Circle		
M-1	C-21	13
J-1	C-22*	17
X-1	C-23*	22
Straight		
SC-1	SC-11*	13
TS	SC-1*	48
KS	SC-2	60
SKS	SC-2	60
SKS-3	SC-6*	90

TROCARPOINT

Ethicon	Medtronic (Covidien)	Needle Length (mm)
Straight		
STP-10	SC-250*	250

Plastic / Cosmetic

PRIME™ REVERSE CUTTING

Ethicon	Medtronic (Covidien)	Needle Length (mm)
1/2 Circle		
P-6	P-16	8
3/8 Circle		
P-1	P-10	11
P-3	P-13	13
PS-3	P-11	16
PS-2	P-12	19
PS-1	P-14	24
PS		26
PSL	P-17	30
PSLX	P-18*	36
1/2 Circle		
P-2	P-21	8
PS-6		10
PS-5	P-22	13
PS-4	P-24	16
PS-4C		16

PC PRIME™ CONVENTIONAL CUTTING

Ethicon	Medtronic (Covidien)	Needle Length (mm)
1/2 Circle		
PC-1	PC13*	13
PC-3	PC-11	16
PC-5	PC-12	19
1/2 Circle		
PC-12		16

PRIME™ REVERSE CUTTING (OPHTHALMIC)

Ethicon	Medtronic (Covidien)	Needle Length (mm)
1/2 Circle		
G-6	HE-6	8
3/8 Circle		
G-1	HE-1	11
G-3	HE-3	13
1/2 Circle		
G-2	HE-2*	8
G-7	HE-7*	8

Needle Comparison Chart

General Surgery

TAPERPOINT

Ethicon	Medtronic (Covidien)	Needle Length (mm)
1/2 Circle		
UCL		13
CT-3	GS-23	22
CT-2	GS-22	26
CT-1	GS-21	36
CT	GS-24	40
CTX	GS-25	48
TP-1	GS-26	65
XLH	GS-27*	70
1/2 Circle - Heavy Body (Mayo)		
MO-7	HGS-23	22
MO-6	HGS-22	26
MO-5	HGS-20*	31
MO-4	HGS-21	36
MO-2	HGS-24	40
1/2 Circle		
UR-6	GU-46	26
UR-5	GU-45	36
UR-4	GU-44	40

Cardiovascular

TAPERPOINT

Ethicon	Medtronic (Covidien)	Needle Length (mm)
1/2 Circle		
BV130-5	MV-135-5/MVF-135-5/ CV-345	6.5
BV175-6	MV175-8/MVF175-8/ CV-351	8
BV175-7		8
BV175-8	MV175-9	9.3
BV-1	CV-1/CVF-1	9.3
BV	CV/CVF/CVL*	11
C-1	CV-11	13
BB	CV-15/CV-327	17
BB-1	CV-13	22
TE		32
1/2 Circle		
TF-6		6.5
TF-4		8
RB-4	CV-20/CV-337	10
TF-1	CV-20/CV-337	10
RB-3	CVF-21*	11
TF	CV-22/CVF-22/CV-330	13
RB-2	CV-22/CVF-22/CV-330	13
RB-1	CV-23/CVF-23/ CV331	17
SH-2	CV-24/CV-316	20
SH-1	CV-2	22
SH	V-20/VF-20/CV-305	26
MH	V-26/CV-300	36
LH		40
1/2 Circle		
TF-5		8.5

- △ Conventional Reverse Cutting Needle
- △ PC PRIME™ Conventional Cutting Needle
- ▽ MICRO-POINT™ Reverse Cutting Needle
- ▽ PRIME™ Reverse Cutting Needle
- ▽ Reverse Cutting Needle
- Bluntpoint Needle
- ETHIGUARD™ Blunted Taperpoint Needle
- △ Conventional Spatula Needle

TAPERCUT™

Ethicon	Medtronic (Covidien)	Needle Length (mm)
1/2 Circle		
V-34	KV-34	36
V-37	KV-37	40
V-40	KV-40	48
V-56	KV-56	60
Straight		
V-47		36
● BLUNTPOINT		
● ETHIGUARD™ BLUNTED TAPERPOINT		
Ethicon	Medtronic (Covidien)	Needle Length (mm)
1/2 Circle		
BP-2		48
BP-1	BP-27	65
BP		80
● ETHIGUARD™ BLUNTED TAPERPOINT		
Ethicon	Medtronic (Covidien)	Needle Length (mm)
1/2 Circle		
SHB	BTV-20	26
CTB-2		26
CTB-1	GBS-21	36
CTB	BGS-24	40
CTXB	BGS-25/BTP-X	48
TPB-1	BGST-29	65

TAPERCUT™

Ethicon	Medtronic (Covidien)	Needle Length (mm)
1/2 Circle		
CC	KV-1/KVF-1/YV-95*	9.3
CC-1	KV-11/KVF-11	13

TAPERCUT™

Ethicon	Medtronic (Covidien)	Needle Length (mm)
1/2 Circle		
V-4	KV-15	17
V-26	KV-26	26
1/2 Circle		
V-5	KV-5	17
V-6		22
V-7	KV-7	26

TAPERPOINT - EVERPOINT™

Ethicon	Medtronic (Covidien)	Needle Length (mm)
1/2 Circle		
BV130-5		6.5
BV175-6		8
BV175-7		8
BV175-8		9.3
BV-1		9.3
BV		11
C-1		13
1/2 Circle		
RB-2		13

TAPERCUT™ - EVERPOINT™

Ethicon	Medtronic (Covidien)	Needle Length (mm)
1/2 Circle		
CC175-8		9.3
CC		9.3
CC-1		13

MICRO-POINT™ Spatula Needle

Center Point Spatula Needle

ULTIMA™ Reverse Spatula Needle

SABRELOC™ Reverse Spatula Needle

TAPERCUT™ Needle

TAPERCUT™ Needle

Taperpoint Needle

Orthopaedic

REVERSE CUTTING

Ethicon	Medtronic (Covidien)	Needle Length (mm)
1/2 Circle		
CP-2	GS-10	26
CP-1	GS-11	36
CP	GS-12	40
CPX	GS-13	48
LS-1		45
1/2 Circle - Heavy Body		
OS-2		17
OS-4	HOS-10	22
OS-6	HOS-11	26
OS-8	HOS-12	40

Endoscopic

TAPERPOINT

Ethicon	Medtronic (Covidien)	Needle Length (mm)
1/2 Circle		
EN-S	SK*	19
3/8 Circle		
EN-2		31

TAPERCUT™

DT-6		20
EN-3		24
EN-4		24
V-38	TJ-32*	30

MicroSurgery

TAPERPOINT - MICRO-POINT™

Ethicon	Medtronic (Covidien)	Needle Length (mm)
1/2 Circle		
BV50-3	MV-50-3	3.8
BV75-3	MV-70-3	3.8
BV75-4	MV-70-4	5.1
BV100-3	MV-100-3	3.8
BV100-4	MV-100-4	5.1
BV130-3	MV-135-3	3.8
BV130-4	MV-135-4	5.1
BV130-5	MV-135-5/MVF-135-5/ CV-345	6.5
1/2 Circle		
BVH100-3		5
Straight		
ST100-3		3

TAPERCUT™ - MICRO-POINT™

Ethicon	Medtronic (Covidien)	Needle Length (mm)
1/2 Circle		
V75-3	MVK-70-3*	3.5
V100-3		3.5
3/8 Circle		
V130-3		3.8
V130-4	MVK-100-4*	4.2
V130-5		5.8

MICRO-POINT™ REVERSE SPATULA

Ethicon	Medtronic (Covidien)	Needle Length (mm)
1/2 Circle		
VAS100-4		5

Sternum

CONVENTIONAL CUTTING

Ethicon	Medtronic (Covidien)	Needle Length (mm)
1/2 Circle		
CCS	SCC	48
CCS-1	SCC-1	60

Ophthalmic

ULTIMA™ REVERSE SPATULA

Ethicon	Medtronic (Covidien)	Needle Length (mm)
1/2 Circle		
CS44-6		6.5
1/4 Circle		
CS90-6		7
CTC-6L		13
3/8 Circle		
CS140-6	SE-140-6	6.5
CS140-8	SE-140-8	6.5
1/2 Circle		
CSM-6		4.3
CS160-6	SE-160-6	5.5
CS175-6	SE-175-6	7
CS160-8	SE-160-8	5.5
CS175-8	SE-175-8	7
Compound Curve		
CSC-6		5

MICRO-POINT™ REVERSE SPATULA

Ethicon	Medtronic (Covidien)	Needle Length (mm)
1/2 Circle		
TG100-8	SE-100-8	6.5
3/8 Circle		
TG140-6	SE-140-6	6.5
TG140-8	SE-140-8	6.5
1/2 Circle		
TG160-4	SE-160-4	5.3
TG160-4-3M		4.3

CONVENTIONAL SPATULA

Ethicon	Medtronic (Covidien)	Needle Length (mm)
1/2 Circle		
SS-29/SE-110-11	S-29	7.6
S-4		8
1/2 Circle		
S-28	SS-28	7.6

SABRELOC™ REVERSE SPATULA

Ethicon	Medtronic (Covidien)	Needle Length (mm)
1/2 Circle		
SM-1	SS-1	6
1/4 Circle		
S-14	SS-14	8
S-24	SS-24	8
1/2 Circle		
S-22	SS-22	8
S-2	S-2	8

SABRELOC™ CENTER POINT SPATULA

Ethicon	Medtronic (Covidien)	Needle Length (mm)
1/2 Circle		
RD-1		8

SABRELOC™ TRANSCAMBER REVERSE SPATULA

Ethicon	Medtronic (Covidien)	Needle Length (mm)
Straight		
STC-6		16

TAPERCUT™

Ethicon	Medtronic (Covidien)	Needle Length (mm)
1/2 Circle		
CIF-4		13

For complete indications, contraindications, warnings, precautions, and adverse reactions, please reference full package insert. The third-party trademarks used herein are trademarks of their respective owners.

*Near equivalent - may have slight variations.
Needle Specialty may vary.
Needles not actual size.
Sales types being transitioned to PRIME™ are categorized under PRIME™ geometry.